

PerCP/Cyanine 5.5 Anti-Mouse/Human Ki-67 antibody

Catalog Number: 202415, 202416
Size: 25 µg, 100 µg

PRODUCT DETAILS

Clone: K67AR2b

Application: Flow Cytometry

Reactivity: Human, Mouse

Format: PerCP/Cyanine5.5

Isotype: Rat IgG2b

Recommended Usage: For flow cytometric staining, it is recommended to use less than 1 µg of this reagent per 0.5-1.0 million cells in a 100 µL volume. Optimal reagent performance should be determined by titration for each specific application. PerCP/Cyanine5.5 has an excitation max at 482 nm and an emission max at 695 nm.